



# **URIT-5380Vet**

Auto 5-Diff Hematology with RET

### All for Vet, All for Better



# **Unique 4-angle Pure Laser Technology**



Without Chemical Dye

Original Cell Morphology Signal

One Channel for WBC 5-Diff

# **URIT-5380Vet**

## Auto 5-Diff Hematology with RET

# **Efficient, Auto-loader**



Up to 50-minute Walk Away

# Reliable, Two 3D Sterograms



- 2 Histograms
- 2 Scatter Grams
- 2 3D Sterogrms

### **Less Maintenance**

Exclusive Software for Veterinary Use

**High Quality Components** 

One-key Maintenance

Integrated Internal Design







### **URIT-5380Vet**

### Auto 5-Diff Hematology with RET

### **Technical Specifications**

#### **Species**

Dog, Cat, Horse, Pig, Cow, Buffalo, Rabbit, Monkey, Rat, Mouse, Sheep, Goat, Camel, 3 User-defined

#### **Parameters**

34 Parameters: WBC, LYM%#, MON%#, NEU%#, EOS%#, BASO%#, RBC, HGB, HCT, MCV, MCH, MCHC, RDW\_CV, RDW\_SD, PLT, MPV, PDW, PCT, P\_LCC, P\_LCR, RET%, RET\_ABS, IRF, NRBC%#, LIC%#, ALY%#

2 Histograms for RBC and PLT, 2 Scattergrams for WBC Diff, 2 Scattergrams for Reticulocyte,

2 3D grams for WBC Diff

### **Principles**

WBC Diff: Laser Light Multi-Dimensional Cell Clasiffication; WBC/HGB/PLT Analysis: Impedance Method; HGB Analysis: Cyanide-free Reagent Colorimetry

#### Reagent

Diluent, Sheath, Detergent, Lyse

#### Sample Position

Up to 50 Samples

#### **Test Mode**

Whole Blood&Prediluted

#### Sample Volume

Whole Blood: 20 µL; Prediluted: 20 µL

### Throughput

Up to 60 T/H

#### **Data Storage**

Up to 200,000 Sample Results with Grams

#### **Display**

PC Display (17 inch Color TFT Screen)

#### **Dimension**

580\*550\*705 mm, 73 kg



Marketed by:

**Uritec India**, Plot No. 2, 1st Floor, NK Sharma Office Road, Lohgarh, Zirakpur -140603, India

Phone: +91 97799 80500 e-mail: info@uritecindia.com Web: www.uritecindia.com